

Output Properties

Properties

General

Name

Output1

Description

Data Type

List of Char

☒ Char

☐ Scalar

☐ Integer

☒ List

☐ Decimal

☐ Boolean

☐ DateTime

☐ Complex

Default Value

Routing Dispositions

User

External

Session

Machine

Remove

Type

Comment

Field	Description
Name	The name of the Output. It can be up to 255 characters long, must begin with a letter or an underscore, and the subsequent characters must be letters, numbers, or underscores. Output names are not case-sensitive.
Description	The description of the Output (localizable).
Data Type	The data type of the Output. The list of valid data types is the following: ▶ Char – character string, ▶ Integer – integer values (to know the range of this Input type, refer to <i>Int32 Struct</i> at Microsoft Docs), ▶ Decimal – decimal values (to know the range of this Input type, refer to <i>Decimal Struct</i> at Microsoft Docs), ▶ Boolean – true or false value, ▶ Date Time – Date and Time (to know the range of this Input type, refer to <i>DateTime Struct</i> at Microsoft Docs), ▶ Complex – an object that can contain any number of properties. Each property can be of a different type (all types are supported, e.g., Char, Integer, List Of). It is also possible that the Complex data type will not have a value (can be treated as null or undefined). You can combine the above types with the ones below: ▶ Scalar – a single data value,

	▶ List – a list of values.
Default Value	The Default Value is used when a Function does not set the value of the Output. This option is disabled when a Composition Group is used.
Routing Dispositions	Outputs may be sent to one or more destinations. The list of valid destinations is as follows: ▶ User ▶ Machine ▶ External Output ▶ Session Variable The Complex/List of Complex data type of the Output can be sent only to Session or External Output destinations. When a Routing Disposition is selected on the list, its properties are displayed below the grid. The Remove Button – removes the selected destination from the field below.